

## cDNA Library, *S. cerevisiae* , Log Phase

#02-701    500 ng

This cDNA library (plasmid DNA) is constructed from *Saccharomyces cerevisiae*, strain S288C- derived poly(A)<sup>+</sup> RNA at the log phase by the Linker-Primer method (Ref.1) by Prof. H. Nojima of Osaka University. This library is unidirectionally cloned by using the oligo (dT)<sub>18</sub> linker primer which contains the restriction enzyme site of *Not*I, and *Bam*HI (*Bgl*II)-*Sma*I adaptor. The pLZ3 vector (shown below) used in this library can not replicate in *S. cerevisiae* but contains pUC ori for replication in *E. coli*

### Application

PCR screening of known or unknown gene: Prepare the primers for the known or unknown gene (cDNA) and amplify the gene by PCR from this library followed by cloning to an appropriate vector. Standard amplifying conditions: 35 cycles of PCR reactions using 10-100 ng of cDNA as a template. (Change the quantity of template and the number of cycles depending on the expression rate of mRNA of the objective gene.)

### Specification

Quantity: 500 ng (40 ng/ul, 13ul) in 10 mM Tris-HCl-1mM EDTA (pH 7.5)

Quality: 1) Number of independent clones:  $3.6 \times 10^6$

2) Average insert size : longer than 1 kb

Storage: -20°C

**References:** Construction of this library is described in Supplementary data of Ref.3

1. Kobori M *et al* "Large scale isolation of osteoclast-specific genes by an improved method involving the preparation of a subtracted cDNA library." *Genes Cells* **3**: 459-475 (1998) PMID: [9753427](#)
2. Tanaka S and Nojima H "Nik1: a Nim1-like protein kinase of *S. cerevisiae* interacts with the Cdc28 complex and regulates cell cycle progression." *Genes Cells* **1**, 905-921 (1996) PMID: [9077450](#)
3. Tougan T, Okuzaki D, Nojima H. Chum-RNA allows preparation of a high-quality cDNA library from a single-cell quantity of mRNA without PCR amplification. *Nucleic. Acids Res.*, 36(15):e92, (2008) PMID:[18603591](#)


### Note

\* This library is to be used only by the purchaser. It is not allowed to amplify and transfer the library to a third person.

\* **Related products:** human tissue specific cDNA libraries and cDNA libraries of model organisms (See HP).

**Distributed by AS ONE International, Inc.**

Fig. Structure of pLZ3 and the restriction sites.



### ; MCS(pLZ3)

```

          CpoI(3)   SauI(b) MluI(5)                AatII(3) BglII(5) AscI(5)       BalI(b)
PstI(3) SacI(3)  ApaI(3)  -----                EcoRI(5) XbaI(5) AflIII(5)       BstXI(5)
SseI(3) -----                T7 Promoter
NNNCTGCA CCTGCAGGAGCTCGGACCGGGCCCTTAGGACGCGTAATACGACTCACTATAGGGAATTCGACGTCTAGATCTTAAGGCGGCCAAGGGGTGGCCA
NNNG  ACGTGGACGTCCTCGAGCCTGGCCCGGAATCCTGCGCATTATGCTGAGTGATATCCCTTAAGCTGCAGATCTAGAATTCGCGCGGTCCCAACCGGT

          BstEII(5)
          -----
          NheI(5)
          -----
SnaBI(b)  DraIII(3) -----  SceI(3)          NotI(5) T3 promoter          SmaI(3)   NruI(b)   SacII(3)
          -----                -----SplI(5)-----  PacI(3)-----  SacI(3)
CGTGGTAACCAACGGGGTGGCTAGCTAGGATAACAGGGTAATATAGCGGCCGCCCTTAGTGAGGGTTAATTTAAATCGTACGTCGGGATTAATTAACCGCGGTGGAGCT CAAT
GCACCAATTGGTGCCCAACCGATCGATCCCTATGTGCCATTTATATCGCGGGGGAAATCACTCCCAATTAATTTAGCATGCAGCGCTAATTAATTGGGCCACC TCGACTTA

TCGCCCTATAGTGAGTCGTATTA -3'
AGCGGGATATCACTCAGCATAAT -5'

```