

## DNA (cytosine-5) methyltransferase (mouse), Dnmt1

#10-201 300 units

DNA methylation is significant for epigenetic regulation of gene expression, X chromosome inactivation, genomic imprinting, and development. Aberrant methylation patterns are associated with certain human tumors and developmental abnormalities. In vertebrates, there are two types of DNA methyltransferase activities; de novo and maintenance types. Two DNA methyltransferases, Dnmt3a and Dnmt3b, are responsible for the creation of methylation patterns at an early stage of embryogenesis (de novo-type), while Dnmt1 is responsible for the maintenance of methylation patterns during replication. Dnmt1 favors to methylate the hemimethylated DNA and preferentially methylates one strand of the double-stranded DNA during its processive methylation. This product, mouse Dnmt1 deleting the N-terminal 290 amino acid residues, was expressed using a baculovirus expression system and purified by Prof. S. Tajima and Dr. I. Suetake of Osaka University (ref.2).

### Uses

- 1) In vitro methylation of cytosine residues in hemimethylated DNA at 5'....CG...3'. (ref. 1, 2)
- 2) Antigen for anti-mammalian Dnmt1 antibodies.

### Specification

Form: 0.5mg protein/ml in 0.2M NaCl, 10mM HEPES (pH 7.4), 50% glycerol

Definition of specific activity: 1 unit is defined as the amount of the enzyme that transfer 1 pmole of methyl group to poly dI-dC substrate during 30 minutes at 37°C

Specific activity: 17 units/ul

Storage: Store at -20°C

Quality Assurance: Greater than 95% protein determined by SDS-PAGE (CBB staining) (Fig.1)

Data Link UniProtKB/Swiss-Prot [P13864](#) (DNMT1\_MOUSE)

---

Distributed by AS ONE International, Inc.

Tel: 408-638-7415

[www.asone-int.com](http://www.asone-int.com)

[info@asone-int.com](mailto:info@asone-int.com)

### Reaction Conditions

Incubate in 1 x Dnmt1 Reaction Buffer (20mM Tris-HCl, pH7.4, 0.5 mM EDTA, 0.2 mM DTT, 5% glycerol) with 10 $\mu$ M S-adenosylmethionine (SAM) at 37°C

### Reagents Supplied with Enzyme

Dnmt1 Reaction Buffer (5 x)

20mM S-adenosylmethionine (SAM) which was purified by chromatography from the commercial reagent and dissolved in H<sub>2</sub>O

Note: SAM is very unstable. Store at -80°C and use it within 6 months.

### References

This product was used in ref.1.

1. Vilkaitis G et al (2005), J Biol Chem 280: 64-72 PMID :

[15509558](https://pubmed.ncbi.nlm.nih.gov/15509558/)

2. Tajima S and Suetake I (1998) , J Biochem 123: 993-999 Review PMID: [9603984](https://pubmed.ncbi.nlm.nih.gov/9603984/)




Fig.1 SDS-polyacrylamide gel electrophoresis of recombinant Dnmt1