

Instruction for the

MSB[®] Spin PCRapace

With only 5 min of processing time, the **MSB[®] Spin PCRapace** is the fastest system for purification of up to 100 µl PCR products (80 bp - 30 kb) from dNTP's, primers, enzymes, additives and salts using the unique **MSB[®] technology** – a washing step is not required.

The kit is useful as well for the cleanup of DNA fragments from salts and enzymes from restrictions digestion, ligation and cDNA synthesis mixtures. For the efficient concentration of PCR products the **MSB[®] Spin PCRapace** can be used as well as for the reliable removal of Dye terminators from DNA cycle sequencing reactions. The recovery of PCR product is 80 – 95 %.

MSB[®] Vario CleanUp

The **MSB[®] Vario CleanUp** provides two different Binding Buffers. Using **Binding Buffer 2** the kit purifies PCR products and removes complete primer dimer simultaneous in a unique way.

Using the **Binding Buffer 1** the kit is a powerful tool for the DNA fragment purifications from diverse enzymatic reaction mixtures, for DNA fragment concentration as well as for the reliable removal of Dye terminators from DNA cycle sequencing reactions. Using the unique **MSB[®] technology** a washing step is not required.

Invisorb[®] Spin DNA Extraction

The **Invisorb[®] Spin DNA Extraction Kit** is the ideal tool for extraction of DNA fragments of 80 bp – 30 kb from standard or low melt agarose gels in TAE and TBE buffers in high end-concentration of DNA. Up to 300 mg agarose gel slices can be processed per spin column.

Invisorb[®] Fragment CleanUp

The **Invisorb[®] Fragment CleanUp** is a combination of different technologies .

Using **MSB[®] technology** the kit is the ideal tool for ultra fast purification of DNA fragments (80 bp – 30 kb) after PCR and other enzymatic reactions, for the efficient concentration of DNA fragments as well as for the reliable removal of Dye terminators from DNA cycle sequencing reactions. The recovery of PCR product is 80 – 95 % and a washing step is not required

Using the **Invisorb[®] technology** the kit is the ideal tool for the efficient extraction of DNA fragments from TAE or TBE agarose gels.

Trademarks: Invisorb[®], Eppendorf[®]. Registered marks, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

The Invisorb[®] technology is covered by patents and patent applications: US 6,110363, US 6,043,354, US 6,037,465, EP 0880535, WO 9728171, WO 9534569, EP 0765335, DE 19506887, DE 10041825.2, WO 0034463.

MSB[®] and Invisorb[®] are registered trademark of STRATEC Molecular GmbH.

The PCR process is covered by US Patents 4,683,195, and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG.

© 2011 STRATEC Molecular, all rights reserved.

Contents

Kit content of the MSB [®] Spin PCRapace	3
Kit content of the MSB [®] Vario CleanUp	3
Kit content of the Invisorb [®] Spin DNA Extraction	4
Kit content of the Invisorb [®] Fragment CleanUp	4
Symbols	5
Storage	5
Quality control	5
Safety information	6
Intended use	7
Product use limitation	7
PCR product & DNA fragment purification and concentration	8
MSB [®] technology	9
Invisorb [®] technology	9
Product characteristics of the MSB [®] Spin PCRapace	10
Product characteristics of the MSB [®] Vario CleanUp	11
Product characteristics of the Invisorb [®] Spin DNA Extraction	12
Product characteristics of the Invisorb [®] Fragment CleanUp	13
Important notes	14
Important points before starting a protocol	14
Equipment and reagents to be supplied by user	14
Principle and procedure of the MSB [®] Kits	15
Scheme MSB [®] Spin PCRapace and MSB [®] Vario CleanUp	16
Principle and procedure of the Invisorb [®] Spin DNA Extraction	17
Scheme Invisorb [®] Spin DNA Extraction	19
Principle and procedure of the Invisorb [®] Fragment CleanUp	20
Scheme Invisorb [®] Fragment CleanUp	21
<i>Protocol 1:</i> Purification and concentration of DNA fragments from enzymatic reactions, like PCR-products from PCR reactions, cDNA synthesis, enzyme restriction digestions	22
<i>Protocol 2:</i> Removal of DyeDeoxy TM terminators from DNA cycle sequencing reactions of PCR-products and Plasmids after use ABI Prism TM terminator Kits	23
<i>Special Protocol 3:</i> Purification and concentration of PCR - products from 200 µl PCR reactions	24
<i>Protocol 4:</i> Purification and concentration of PCR-products over 100 bp with contaminations of primer dimers or shorter fragments*	25
<i>Protocol 5 :</i> Extraction of a DNA-fragment from an agarose gel slice	26
<i>Protocol 6 :</i> Instruction for repair of A-overhangs in DNA fragments after gel purification	26
Troubleshooting	27
Ordering information	28

Kit contents of the MSB[®] Spin PCRapace

Store all kit components at room temperature

	10 preps	50 preps	250 preps	500 preps
Catalogue No	1020220900	1020220200	1020220300	1020220400
Binding Buffer	6 ml	30 ml	140 ml	2 x 140 ml
Elution Buffer	2 ml	3 x 2 ml	30 ml	60 ml
Spin Filter	10	50	250	500
2.0 ml Receiver Tubes	10	50	250	500
1.5 ml Receiver Tubes	10	50	250	500
Manual	1	1	1	1

Kit contents of the MSB[®] Vario CleanUp

Store all kit components at room temperature

	10 preps	50 preps	250 preps	500 preps
Catalogue No	1020230900	1020230200	1020230300	10202230400
Binding Buffer 1	6 ml	30 ml	140 ml	2 x 140 ml
Binding Buffer 2*	6 ml	30 ml	140 ml	2 x 140 ml
Elution Buffer	2 ml	3 x 2 ml	30 ml	60 ml
Spin Filter	10	50	250	500
2.0 ml Receiver-Tubes	10	50	250	500
1.5 ml Receiver-Tubes	10	50	250	500
Manual	1	1	1	1

* Binding Buffer 2 is only stable for 6 months!

Kit contents of the Invisorb® Spin DNA Extraction Kit

Store all kit components at room temperature

	3 preps	10 preps	50 preps	250 preps
Catalogue No	1020110100	1020110900	1020110200	1020110300
Gel Solubilizer S	2 x 2 ml	12 ml	60 ml	2 x 140 ml
Binding Enhancer	2 ml	6 ml	30 ml	3 x 50 ml
Wash Buffer	15 ml (ready to use)	15 ml (ready to use)	18 ml (final volume 60 ml)	2 x 45 ml (final volume 2 x 150 ml)
Elution Buffer	2 ml	2 ml	2 x 2 ml	15 ml
Spin Filter	3	10	50	250
2.0 ml Receiver Tubes	3	10	50	250
1.5 ml Receiver Tubes	3	10	50	250
Manual	1	1	1	1
Initial steps			add 42 ml 96-100% ethanol to the bottle Wash Buffer	add 105 ml 96-100% ethanol to each bottle Wash Buffer









Kit contents of Invisorb® Fragment CleanUp

Store all kit components at room temperature

	10 preps	50 preps	250 preps
Catalog No	1020300900	1020300200	1020300300
Gel Solubilizer S	12 ml	60 ml	2 x 140 ml
Binding Buffer	6 ml	30 ml	140 ml
Binding Enhancer	6 ml	30 ml	150 ml
Wash Buffer	15 ml (ready to use)	18 ml (final volume 60 ml)	2 x 45 ml (final volume 2 x 150 ml)
Elution Buffer	2 ml	2 x 2 ml	15 ml
Spin Filter	10	50	250
2.0 ml Receiver Tubes	10	50	250
1.5 ml Receiver Tubes	10	50	250
Manual	1	1	1

Initial steps		add 42 ml 96-100% ethanol to the bottle Wash Buffer	add 105 ml 96-100% ethanol to each bottle Wash Buffer
----------------------	--	---	--

Symbols

	manufacturer
	lot number
	catalogue number
	date of manufacture
	expiry date
	consult operating instructions
	temperature limitation
	do not reuse

Storage

All buffers of the **MSB[®] Spin PCRapace**, **MSB[®] Vario CleanUp**, **Invisorb[®] Spin DNA Extraction Kit** and the **Invisorb[®] Fragment CleanUp** should be stored well sealed and dry at room temperature and are stable for at least 12 months under these conditions, except **Binding Buffer 2** (only stable for 6 months!).

Before every use make sure that all components are at room temperature. If there are any precipitates within the provided solutions dissolve these precipitates by warming up carefully.

Room temperature (RT) is defined as range from 15 - 30°C.

Quality control

STRATEC Molecular guarantees the correct function of the **MSB[®] Spin PCRapace**, **MSB[®] Vario CleanUp**, **Invisorb[®] Spin DNA Extraction Kit** and the **Invisorb[®] Fragment CleanUp** for applications as described in the manual.

In accordance with STRATEC Molecular's certified Quality Management System all components of the **MSB[®] Spin PCRapace**, **MSB[®] Vario CleanUp**, **Invisorb[®] Spin DNA Extraction Kit** and the **Invisorb[®] Fragment CleanUp** were tested against predetermined specifications to ensure consistent product quality.

If you have any questions or problems regarding any aspects of the **MSB[®] Spin PCRapace**, **MSB[®] Vario CleanUp**, **Invisorb[®] Spin DNA Extraction Kit** and the **Invisorb[®] Fragment CleanUp** or other STRATEC Molecular products, please do not hesitate to contact us.

For technical support or further information please contact:

from Germany +49-(0)30-9489-2901/ 2910

from abroad +49-(0)30-9489-2907

or contact your local distributor.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. Avoid skin contact! Always follow the legal requirements for working with biological material!

For more information, please consult the appropriate material safety data sheets (MSDS). These are available online in convenient and compact PDF format at www.invitek.de for each STRATEC Molecular kit and kit component.

If the buffer bottles are damaged or leaking, wear gloves, and protective goggles when discarding the bottles in order to avoid any injuries.

STRATEC Molecular has not tested the liquid waste generated by the **MSB[®] Spin PCRapace**, **MSB[®] Vario CleanUp**, **Invisorb[®] Spin DNA Extraction Kit** and the **Invisorb[®] Fragment CleanUp** procedures for residual infectious materials. Contamination of the liquid waste with residual infectious materials is highly unlikely, but cannot be excluded completely. Therefore, liquid waste must be considered infectious and must be handled and discarded according to local safety regulation.

Below are listed the European Community risk and safety phrases for the components of the **MSB[®] Spin PCRapace**, **MSB[®] Vario CleanUp**, **Invisorb[®] Spin DNA Extraction Kit** and the **Invisorb[®] Fragment CleanUp** to which they apply.

Binding Buffer 1



danger
H225-319-336 P210-233-305-351-338

Binding Buffer 2



danger
H225-319-336 P210-233-305-351-338

Binding Enhancer



danger
H225-319-336 P210-233-305-351-338

Gel Solubilizer



warning
H302-312-332-412 EUH032 P273

H225:	Highly flammable liquid and vapour.
H319:	Causes serious eye irritation.
H336:	May cause drowsiness or dizziness.
H302:	Harmful if swallowed.
H312:	Harmful in contact with skin.
H332:	Harmful if inhaled.
H412:	Harmful to aquatic life with long lasting effects.
EUH032:	Contact with acids liberates very toxic gas.
P210:	Keep away from heat/sparks/open flames/hot surfaces. — No smoking.
P233:	Keep container tightly closed.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P273:	Avoid release to the environment.

Emergency medical information in english and german language can be obtained 24 hours a day from:

Poison Information Center Freiburg, Germany: Tel.: ++49 761-19240

Intended use

See for each product the product characteristic, there are presented the intended use and the product use limitation for each kit!

These products are intended for use by professional users such as technicians, physicians and biologists trained in molecular biological techniques. They are designed to be used with any downstream application employing enzymatic amplification or other enzymatic modification of DNA followed by signal detection or amplification. Any diagnostic results generated using the sample preparation procedure in conjunction with any downstream diagnostic assay should be interpreted with regard to other clinical or laboratory findings.

To minimize irregularities in diagnostic results, adequate controls for downstream applications should be used.

The kits are developed, designed, and sold for research purposes only.

Product use limitation

For purification the fragments should not contain less than 70 bp,

The maximum length of primers which can be removed is 40 bp.

DNA fragments should not be bigger than 30 kb and not smaller than 70 bp.

The included chemicals are only useable once.

There is no guarantee of functionality if using non-recommended starting material or flow trace.

The user is responsible for validating the performance of the STRATEC Molecular kits for any particular use, since the performance characteristics of our kits have not been validated for any specific application. STRATEC Molecular kits may be used in clinical diagnostic laboratory systems after the laboratory has validated the complete diagnostic system as required by CLIA' 88 regulations in the U.S. or equivalents in other countries.

All products sold by STRATEC Molecular are subjected to extensive quality control procedures (according to ISO 9001-2000 and ISO EN 13485) and are warranted to perform as described when used correctly. Any problems should be reported immediately.

The chemicals and the plastic parts are for laboratory use only, they must be stored in the laboratory and must not be used for purposes other than intended.

The kit contents are unfit for consumption.

PCR product & DNA fragment purification and concentration

This manual characterizes innovative kits for membrane adsorption based purification of PCR products or DNA fragments from enzymatic reaction mixtures and of DNA fragments from agarose gels using the high-performance **MSB®** or the well established **Invisorb® technology**. Different kits were developed and adapted to suit different DNA fragment purification needs.

Sample Source	Recommended Kit***	Sample Volume	Recovery
PCR products, ligation reaction mixture, enzyme digestion mixture, cDNA synthesis mixture cycle sequencing reaction mixture	MSB® Spin PCRapace Invisorb® Fragment CleanUp MSB® Vario Cleanup	up to 100 µl	80 – 95 %*
PCR products with primer dimer cycle sequencing reaction mixture with dye terminators	MSB® Vario Cleanup	up to 100 µl	80 – 95 %*

Sample Source	Recommended Kit***	Sample Volume	Recovery
Strong contaminated DNA, DNA fragments and PCR products	Invisorb® DNA CleanUp (see ordering information page 28)	up to 200 µl	60 - 85 %*
Agarose gels (TAE, TBE)	Invisorb® Spin DNA Extraction Invisorb® Fragment CleanUp	up to 300 mg gel slices	60 – 90 %**

Advantages:

- convenient and fast sample processing
- the most efficient removal of contaminants
- high recovery rate of PCR products or DNA fragments
- broad range of fragment sizes 80 bp - 30 kb can be purified

The purification procedure has been optimized to remove salts, enzymes, nucleotides, mineral oil, agarose, ethidium bromide and other impurities from DNA samples. Specialized binding buffers promote selective adsorption of DNA fragments and PCR products. The pure DNA is eluted in a small volume of buffer or water, ready to use for any subsequent application. The innovative **MSB® technology** is characterized on the following pages.

MSB® technology

represents the fastest technology for purification of DNA fragments with high recovery rates.

This development opens new possibilities for simplifying the purification procedure and to considerably reduce handling steps and processing time. The **MSB® kits** have been designed for extremely efficient purification and/or concentration of PCR products or of DNA fragments from enzymatic reaction mixtures with **only two steps**.

The DNA fragments bind to the membrane at minimal concentrations of non chaotropic salts. Therefore a washing step is not required. High concentrated, pure DNA fragments are eluted ready to use.

Advantages:

- ultra fast and easy (two step format), only binding and elution
- excellent purity without washing
- 80 – 95 % rate of recovery

Invisorb® technology

is the first technology for the extraction of highly purified nucleic acids using non chaotropic binding conditions.

Starting from complex biological samples, the method allows binding of nucleic acids to nearly all surfaces, such as membranes, carrier or magnetic particles. In combination with specially modified surfaces a selective binding of different nucleic acid targets can be realized. Binding of nucleic acids under non chaotropic salt conditions leads to

- high binding efficiency
- reproducible high yields of ready to use DNA simplified protocols; reliable, time saving and easy handling steps
- guaranteed freedom of operation for all fields of applications
- improved quality assurance in the customers laboratories

Product characteristics of the MSB[®] Spin PCRapace

Starting material	Yield	Time for preparation
up to 100 µl reaction volume like PCR reaction mixture, up to 100 µl restriction digestion mixture, up to 100 µl ligation mixture, up to 100 µl cDNA synthesis mixture, up to 100 µl cycle sequencing reaction mixture	80 – 95 %, depends on fragment length	approx. 5 min

The **MSB[®] Spin PCRapace** provides a convenient tool for ultra fast and efficient direct purification of about 100 µl PCR products from 80 bp up to 30 kb from amplification reactions. With max. 5 min of processing time the **MSB[®] Spin PCRapace** is the fastest system for the separation of PCR products (80 bp - 30 kb) from dNTP's, primers, enzymes, additives, labelling reagents (biotin, radioactive ATP etc.) and salts. The recovery of PCR product is 80 – 95 %. The kit is further useful for DNA fragment cleanup from

- restriction digestion mixture
- ligation mixture
- cDNA synthesis mixtures

The kit is also a powerful and efficient tool for the

- concentration of DNA fragments as well as
- for purification of linearized pDNA from restriction mixtures and
- for the reliable removal of Dye terminators from DNA cycle sequencing reactions

The DNA-fragments will be bound directly onto the surface of a spin filter column based on new buffer composition. No additional and common used washing steps are necessary. Finally, the DNA fragments will be eluted with low salt buffer or H₂O.

MSB [®] Spin PCRapace	parameters
Binding capacity for fragment size of 200 – 500 bp Binding capacity for fragment size of 20.000 kb Elution volume	up to 10 µg up to 100 µg minimum 10 µl

Beside the extremely time efficient procedure the kits provides high yield and purity of the recovered DNA fragments. The recovery rate (80 - 95 %) is nearly independent from PCR additives, the PCR reaction volume and the fragment size. The purified PCR product can be used in subsequent downstream applications:

- sequencing
- cloning and labeling experiments
- hybridization
- transcription
- digestion with restriction enzymes

Product use limitation

The kit works not suitable with PCR products smaller than 70 bp, single stranded DNA may diverge.

Product characteristics of the MSB[®] Vario CleanUp

Starting material	Yield	Time for preparation
up to 100 µl reaction volume like PCR reaction mixture, up to 100 µl reaction mixtures with primer dimer contamination up to 100 µl restriction digestion mixture, up to 100 µl ligation mixture, up to 100 µl cDNA synthesis mixture, up to 100 µl cycle sequencing reaction mixture	80 – 95 %, depends on fragment length	approx. 5 min

The speciality of the **MSB[®] Vario CleanUp** is that it provides two different **Binding Buffers**. By using the Binding Buffer 2 the kit:

- PCR product
- removes complete **primer dimer** in a unique way.
- ligation rate of primer dimer free PCR products is > 90%

The **MSB[®] Vario Cleanup** provides further a convenient tool for ultra fast and efficient direct purification of

- PCR products from 80 bp up to 30 kb from amplification reactions,
- DNA fragments from diverse enzymatic reaction mixtures
- for reliable removal of DyeDeoxy[™] terminators from DNA-cycle sequencing reactions (Binding Buffer 1).

Using the unique **MSB[®] technology** a washing step is not required. The kit is also a powerful and efficient tool for the

- concentration of DNA fragments as well as

The DNA-fragments will be bound directly onto the surface of a spin filter column based on new buffer composition. No additional and common used washing steps are necessary. Finally, the DNA fragments will be eluted with low salt buffer or H₂O.

MSB [®] Vario CleanUp	parameters
Binding capacity for fragment size of 200 – 500 bp Binding capacity for fragment size of 20.000 kb Elution volume	up to 10 µg up to 100 µg minimum 10 µl

The isolation protocol as well as all buffers are optimized to provide high yield and purity of the recovered DNA fragment and for effective removal of free and labelled dNTP's from fluorescent reaction mixtures. The primer dimer free purified PCR product can be used in subsequent downstream applications:

- cloning and labeling experiments
- ligation
- sequencing
- hybridization
- transcription
- digestion with restriction enzymes

Product use limitation

The kit works not suitable with PCR products smaller than 70 bp , single stranded DNA may diverge.

Product characteristics of the Invisorb® Spin DNA Extraction Kit

Starting material Size of DNA fragments: 80 bp - 30 kb	Yield	Time for preparation
up to 300 mg of gel slices of from TAE or TBE buffered systems	70 - 90 %, depends on fragment length and kind of agarose gel	approx. 20 min

The **Invisorb® Spin DNA Extraction Kit** provides a convenient tool for fast and efficient extraction and purification of DNA fragments directly from TAE and TBE agarose gels slices in less than 20 min in the comfortable spin filter format. The purification procedure removes agarose, ethidium bromide and other impurities from the DNA sample.

After gel solubilization DNA fragments from 80 bp up to 30 kb are bound directly onto the membrane of a spin filter column. All impurities are removed very efficiently by a washing step. The DNA is eluted in as little as 20 µl low-salt buffer.

Invisorb® Spin DNA Extraction Kit	parameters
Binding capacity for fragment size of 200 – 500 bp Binding capacity for fragment size of 20.000 kb Elution volume	up to 10 µg up to 100 µg minimum 20 µl

The handling and the composition of the **Invisorb® Spin DNA Extraction Kit** is optimized to guarantee high recovery and purity of the DNA fragments. The purified DNA fragments are ready to use for various downstream applications

- digestion with restriction enzymes
- hybridization, labeling, cloning
- sequencing
- *in vitro* Transcription

Product use limitation

It is not suitable for purification of circular plasmids, because as a result of their very different configurations they were detected in the gel in many positions. Furthermore purification of linearized pDNA is possible.

Important Notes

If processing PCR with Taq-Polymerase, possibly the “a overhangs” can become lost during the extraction. It is recommended to perform a reparation of a overhangs subsequent to the extraction if they are needed in the following downstream application. (see page 26)

Product characteristics of the Invisorb[®] Fragment CleanUp

The Invisorb[®] Fragment CleanUp is a combination of MSB[®] Spin PCRapace and Invisorb[®] Spin DNA Extraction Kit.

Starting material Size of DNA fragments: 80 bp - 30 kb	Yield	Time for preparation
up to 100 µl reaction volume like PCR reaction mixture, up to 100 µl restriction digestion mixture, up to 100 µl ligation mixture, up to 100 µl cDNA synthesis mixture, up to 100 µl cycle sequencing reaction mixture	80 – 95 %, depends on fragment length	approx. 5 min
up to 300 mg of gel slices of from TAE or TBE agarose gels	70 - 90 %, depends on fragment length and kind of agarose gel	approx. 20 min

The Invisorb[®] Fragment CleanUp provides ultra fast purification and concentration of up to 100 µl PCR-products from 80 bp up to 30 kb and of other enzymatic reaction mixtures, or of linearized pDNA as well as for the purification of DNA-fragments from agarose gels.

For PCR cleanup the MSB[®] procedure is offered -no commonly used washing steps are necessary. Finally, the DNA fragments will be eluted with low salt buffer or H₂O.

For purification of DNA-fragments from agarose gels the Invisorb[®] procedure is offered. The DNA fragments are bound directly onto the surface of a spin filter column after gel solubilization. The DNA – fragments will be eluted in a low salt buffer after washing.

The extraction protocol as well as all buffers are optimized to provide high yield and purity of the recovered DNA-fragment. The “hands-on time” necessary for the whole procedure is reduced to a minimum. The purification process will be ready in 5 - 20 minutes. The purified DNA-fragments are ready to use in various downstream applications such as:

- digestion with restriction enzymes
- hybridization, labelling, cloning
- sequencing
- *In vitro* Transcription
- ligation and transformation
- DNA sequencing
- amplification, microinjection

Product use limitation

It is not suitable for extraction of circular plasmids, because as a result of their very different configurations they were detected in the gel in many positions.

If processing PCR with Taq-Polymerase, possibly the “a overhangs” can become lost during the extraction. It is recommended to perform a reparation of a overhangs subsequent to the extraction if they are needed in the following downstream application. (see page 26)

Important notes

Important points before starting a protocol

After receiving the kit, check the kit components for damage. If buffer bottles are damaged, contact the STRATEC Molecular Technical Services or your local distributor. In case of liquid spillage, refer to "Safety information" (page 6). Do not use damaged kit components, since their use may lead to poor kit performance.

- always change pipet tips between liquid transfer. To avoid cross-contamination, recommend the use of aerosol-barrier pipet tips
- all centrifugation steps are carried out at room temperature
- when working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles
- discard gloves if they become contaminated
- do not mix kit components with components from other kits unless the lot numbers are identical
- avoid microbial contamination of the kit reagents
- to minimize the risk of infections from potentially infectious material, we recommend working under laminar air-flow until the samples are lysed
- this kit should only be used by personnel trained in vitro diagnostic laboratory practice

Equipment and reagents to be supplied by user

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDS). (See our webpage : www.invitek.de)

- micro centrifuge ($\geq 13.400 \times g$ (12.000 rpm))
- ethanol (96-100%)
- thermo shaker
- pipettes and filter tips
- scalpel 120
- 1.5 ml and 2.0 ml reaction tubes

Principle and procedure of the MSB[®] Kits

The MSB[®] Spin PCRapace & MSB[®] Vario CleanUp procedure comprises the following steps:

- selective binding of DNA fragments to the surface of the DNA-Binding Spin Filter
- elimination of contaminants and ethanol.
- elution of highly pure DNA fragment or PCR product

Binding of DNA fragments

The reaction mixture is mixed in a ratio of 1:5 with the **Binding Buffer** to provide the appropriate condition for the binding of DNA fragments in range of 70 bp - 30 kb to the silica membrane under minimal concentrations of non chaotropic salts.

The binding of small DNA fragments can be supported by the addition of small amounts of isopropanol, but this ratio is very sensitive.

Removal of Contaminants

The DNA fragments bind to the membrane at minimal concentrations of non chaotropic salts. Therefore a washing step is not required. Unwanted primers and impurities, such as salts, enzymes, unincorporated nucleotides, dyes, ethidium bromide, oils, and detergents do not bind to the silica membrane; instead, are drawn through the column by centrifugal force together with the big excess of Binding Buffer. Any remaining Binding Buffer, which may interfere with subsequent enzymatic reactions, is removed by additional centrifugation.


Elution of PCR products or DNA fragments

The efficiency of elution strongly depends on the concentration of salt and on the pH of the elution buffer. Unlike absorption, elution is most efficient under basic conditions and low concentrations of salt. High concentrated, pure DNA fragments are eluted with a minimum of 10 µl of Tris buffer or water. Maximum elution efficiency is achieved between pH 7.0 and 8.5; therefore, when using water to elute, be certain that the pH is within this range. DNA must also be stored at -20 °C when eluted with water, as it is possible that DNA may degrade without a buffering agent.

Sampling and storage of starting material

Best results are obtained using freshly prepared PCR or enzymatic reaction mixtures to prevent DNA digestion. The samples can be stored for some weeks at 4°C .

Scheme of MSB[®] Spin PCRapace & MSB[®] Vario CleanUp



Principle and procedure of the Invisorb® Spin DNA Extraction Kit

The Invisorb® Spin DNA Extraction Kit procedure comprises the following steps:

- excise of DNA-fragments from the agarose gel with a sharp scalpel
- gel removal and binding of DNA fragments on the membrane of the spin column
- selective binding of DNA fragments to the surface of the DNA-Binding Spin Filter
- elimination of contaminants and ethanol
- elution of highly pure DNA fragment

Sampling and storage of starting material

Best results are obtained using freshly prepared DNA slices from gel. Do not expose the gel with the DNA fragment to UV light for a long time. Reduce cutting time under UV light to a minimum. It is damaging to the DNA. The samples can be stored for some weeks at 4°C .

Excise of DNA fragments

For best results the pieces of gel should be small as possible. Use low melting and standard gels with TAE or TBE buffer systems.

Gel removal and binding of DNA fragments

Gel Solubilizer S in the Invisorb® Spin DNA Extraction Kit solubilizes the agarose gel slice under higher temperatures. Together with Binding Enhancer and Binding buffer, it provides the appropriate condition for the binding of DNA to the silica membrane under high concentrations of salt.

Removal of contaminants and of ethanol

The DNA fragments bind to the membrane contaminants and salts are washed away by the ethanol-containing Wash Buffer. Any remaining Wash Buffer which may interfere with subsequent enzymatic reactions, is removed by additional centrifugation.

Elution of DNA fragments

The efficiency of elution strongly depends on the concentration of salt and on the pH of the elution buffer. Unlike to absorption, elution is most efficient under basic conditions and low concentrations of salt. High concentrated, pure DNA fragments are eluted with a minimum of 10 µl of Tris buffer or water. Maximum elution efficiency is achieved between pH 7.0 and 8.5; therefore, when using water to elute, be certain that the pH is within this range. DNA must also be stored at -20 °C when eluted with water, as it is possible that DNA may degrade without a buffering agent.

DNA yield and concentration

DNA yield is dependent on the following three factors: the volume of elution buffer, how the buffer is applied to the column, and the incubation time of the buffer on the column. Between 100 and 200 µl of the elution buffer completely covers the Invisorb® membrane, ensuring maximum yield, even when not applied directly to the center of the membrane. Elution with up to 50 µl requires the buffer to be added directly to the center of the membrane, and an additional 1-2 minute incubation is required for maximum yield if elution is performed with the minimum recommended volume of 10 µl.

DNA fragment purification from agarose gel

- Gel removal and binding of DNA fragments on the membrane of the spin column
- Removal of contaminants and elimination of ethanol
- Elution of DNA fragments

Gel removal and binding of DNA fragments

Gel Solubilizer S in the **Invisorb® Fragment CleanUp Kit** solubilizes the agarose gel slice and, together with the binding enhancer and binding buffer, provides the appropriate condition for the binding of DNA to the silica membrane under high concentrations of salt.

Removal of Contaminants

The DNA fragments bind to the membrane at minimal concentrations of non chaotropic salts. Therefore a washing step is not required. Unwanted primers and impurities, such as salts, enzymes, unincorporated nucleotides, agarose, dyes, ethidium bromide, oils, and detergents do not bind to the silica membrane; instead, they flow through the column. Salts are washed away by the ethanol-containing Wash-Buffer PE. Any remaining Wash Buffer, which may interfere with subsequent enzymatic reactions, is removed by additional centrifugation.


Elution of DNA fragments

The efficiency of elution strongly depends on the concentration of salt and on the pH of the elution buffer. Unlike absorption, elution is most efficient under basic conditions and low concentrations of salt. High concentrated, pure DNA fragments are eluted with a minimum of 10 µl of Tris buffer or water. Maximum elution efficiency is achieved between pH 7.0 and 8.5; therefore, when using water to elute, be certain that the pH is within this range. DNA must also be stored at –20 °C when eluted with water, as it is possible that DNA may degrade without a buffering agent.

DNA yield and concentration

DNA yield depends on the following three factors: the volume of elution buffer, how the buffer is applied to the column, and the incubation time of the buffer on the column. Between 100 and 200 µl of the elution buffer completely covers the **Invisorb®** membrane, ensuring maximum yield, even when not applied directly to the center of the membrane. Elution with up to 50 µl requires the buffer to be added directly to the center of the membrane, and an additional 1-2 minute incubation is required for maximum yield if elution is performed with the minimum recommended volume of 10µl.

General Scheme of Invisorb® Spin DNA Extraction Kit

	<p>Please read the protocols carefully prior to the start of the preparation procedure</p> <hr/> <p>Transfer gel slices (max. 300 mg) into a 1.5 ml or 2.0 ml microcentrifuge tube (not provided)</p> <p>to gel slices up to 150 mg add 500 µl Gel Solubilizer S to gel slices > 150 mg add 1 ml of Gel Solubilizer S and incubate at 50°C for 10 minutes until the gel is completely solubilized</p> <p>add 250 µl Binding Enhancer to a 500 µl reaction volume add 500 µl Binding Enhancer to a 1 ml reaction volume and mix by pipetting up and down for 2-3 times</p> <p>load app. 800 µl of the sample onto the Spin Filter. centrifuge at 10.000 – 12.000 x g (10.500 – 11.300 rpm) for 1 min discard the filtrate for reaction volumes > 800 µl reload the residual volume and repeat the centrifugation step</p> <p>add 500 µl Wash Buffer to the Spin Filter centrifuge for 30 sec at 8.000 x g (9.000 rpm) discard the filtrate. repeat the washing step once again centrifuge for 30 sec at 8.000 x g (9.000 rpm) discard the filtrate</p> <p>remove the residual ethanol by centrifugation for 4 min at maximum speed</p> <p>transfer the Spin Filter into a new 1.5 ml Receiver Tube. add at least 20 µl Elution Buffer directly onto the center of the Spin Filter incubate at room temperature for 5 min centrifuge for 1 minute at 8.000 x g (9.000 rpm) DNA is now ready to use</p>
--	---

Principle and procedure of the Invisorb® Fragment CleanUp

The Invisorb® Fragment CleanUp Kit combines the convenience of spin-column technology with the selective bonding properties of a uniquely designed silica membrane. The Invisorb® columns are designed to give high end-concentrations of purified DNA fragments for subsequent reactions. Special buffers provided with each kit are optimized for efficient recovery of DNA and removal of contaminants in each specific application.

- DNA fragment purification from agarose gel
- Cleanup of PCR mixture or of enzymatic reaction mixtures
- Two handling options: spin columns can be processed in a micro centrifuge or on an Invisorb® 96 Vacuum Manifold

During DNA cleanup using the **Invisorb® Fragment CleanUp Kit**, all enzymes are removed, independent of size and secondary structure.

Protein	Molecular weight per enzyme subunit (kDa)
DNA Polymerase I	109
Klenow fragment	62
Calf intestinal alkaline phosphatase (CIP)	69
T4 DNA ligase	55
T4 Polynucleotide kinase	35
Terminal transferase	32
DNase I	31

No additional and commonly used washing and drying steps are necessary. Finally, the DNA fragments will be eluted with low salt buffer or H₂O.


For purification of DNA-fragments from agarose gels, the DNA fragments from 80 bp up to 30 kb will be bound directly onto the surface of a spin filter column after gel solubilization. The DNA – fragments will be eluted in a low salt buffer after washing.

The extraction protocol as well as all buffers are optimized to provide high yield and purity of the recovered DNA-fragment. The “hands-on time“ necessary for the whole procedure is reduced to a minimum. The purification process will be ready in 5 - 20 minutes.

The purified DNA-fragments are ready to use in various downstream applications such as:

- Digestion with restriction enzymes
- Hybridization
- Labelling
- Cloning
- Sequencing
- *In vitro* Transcription
- Ligation and transformation
- Radioactive and fluorescent sequencing
- Amplification
- Microinjection

General Scheme of Invisorb® Fragment CleanUp



Protocol 1: Purification and concentration of DNA fragments from enzymatic reactions , e.g. PCR-products from PCR reactions, cDNA synthesis, enzyme restriction digestions

Please read the instructions carefully and conduct the prepared procedure !

for MSB[®] Spin PCRapace, MSB[®] Vario CleanUp and Invisorb[®] Fragment CleanUp

Note: *Before starting with the purification procedure please place a Spin Filter into a 2.0 ml Receiver Tube !*

1. Binding of the PCR or DNA - fragments

A. For PCR-mixtures up to 50 µl

Add **250 µl Binding Buffer/ Binding Buffer 1** to the PCR sample and mix very well by pipetting or vortexing. Transfer the sample completely onto a Spin Filter and centrifuge for 3 min at 13.400 x g (12.000 rpm).

B. For PCR-mixture > 50 µl up to 100 µl

Add **500 µl Binding Buffer/ Binding Buffer 1** to the PCR sample and mix very well by pipetting or vortexing. Transfer the sample completely onto a Spin Filter and centrifuge for 1 min at 13.400 x g (12.000 rpm). Remove the filtrate and centrifuge again for 2 minutes.

2. Elution of the PCR or DNA - fragments

Place the Spin Filter into a new 1.5 ml Receiver Tube.

Add at least 10 µl Elution Buffer (or ddH₂O) directly onto the center of the Spin Filter.

Incubate for 1 minute at room temperature. Centrifugation for 1 minute at 9.300 x g (10.000 rpm).

Important Notes:

1. If the PCR-mixture contains mineral oil, we recommend the addition of 500 µl of Binding Buffer independent of the starting volume. It is also possible to wash the bound PCR-fragment once with 500 µl of Binding Buffer.*

2. To increase the final DNA yield we recommend an extended incubation time with Elution Buffer (up to 5 minutes), which will lead to a slightly higher final yield.

3. For concentration of PCR-fragments it is possible to elute with lower volume of Elution Buffer, than the volume of the starting PCR-mixture. The minimum volume is 10 µl.

4. For ligation mixtures please note, that ligation reactions give very often non wanted side products. These also are purified and enriched.

*** You can order additional Binding Buffer (Cat. No. 10202220) for this application by phone (+49 (0)30 948928910/ 01).**

Protocol 2: Removal of DyeDeoxy™ terminators from DNA cycle sequencing reactions of PCR-products and Plasmids after use ABI Prism™ terminator Kits

Please read the instructions carefully and conduct the prepared procedure !

for MSB® Spin PCRapace, MSB® Vario CleanUp and Invisorb® Fragment CleanUp

Note: *Before starting with the purification procedure please place a Spin Filter into a 2.0 ml Receiver Tube !*

1. Binding of the (fluorescent) labeled DNA

Add **500 µl Binding Buffer** to the completed cycle sequencing reaction (20 – 100 µl) and mix very well by pipetting or vortexing. Transfer the sample completely onto a Spin Filter and centrifuge for 4 min at maximum speed.

Attention: *Optional Step 1 if you need primer close by sequence, than perform step one like described below*

This step may help if sequence has to be readable close up to the starting Oligonucleotide but may reduce purity in some reactions (Dye Blubs may appear). Also sequence reactions in this region, depending on the sequence, sometimes show a bad performance.

Add 500 µl **Binding Buffer** to the completed cycle sequencing reaction (20 – 100 µl). Add 150 µl of **Isopropanol** to the mixture and mix very well by pipetting or vortexing. Transfer the sample completely onto a Spin Filter and centrifuge for 4 min at maximum speed. In case of removal of blobs coming from contamination by primer dimers the addition of Isopropanol is not helpful and should be omitted.

2. Elution of the (fluorescent) labeled DNA

Place the Spin Filter into a new 1.5 ml Receiver Tube. Add at least 10 µl Elution Buffer (or ddH₂O) directly onto the center of the Spin Filter. Incubate for 1 minute at room temperature. Centrifuge for 1 minute at 10.000 rpm. Discard the Spin Filter and proceed with the ABI sample loading.

Special Protocol 3: Purification and concentration of PCR - products from 200 µl PCR reactions

Please read the instructions carefully and conduct the prepared procedure !

for MSB® Spin PCRapace, MSB® Vario CleanUp and Invisorb® Fragment CleanUp

Note: *Before starting with the purification procedure please place a Spin Filter into a 2.0 ml Receiver Tube !*

1. Binding of the PCR-fragments

For PCR-mixture 200 µl

Add **1000 µl Binding Buffer** to the PCR sample and mix very well by pipetting or vortexing. Transfer the sample in two aliquots onto a Spin Filter and centrifuge for 3 min at 12.000 rpm. Remove the filtrate and centrifuge again for 2 minutes

2. Elution of the PCR-fragments

Place the Spin Filter into a new 1.5 ml Receiver Tube.
Add at least 10 µl Elution Buffer (or ddH₂O) directly onto the center of the Spin Filter.

Important Notes:

- 1. The provided volume of Binding Buffer is calculated based on the needed buffer volumes in protocol 1 and 2. The needed amount of protocol 3 is not considered.*
- * You can order additional Binding Buffer (Cat. No. 10202220) for this application by phone (+49 (0)30 948928910/ 01) or by mail.*
- 2. If the PCR-mixture contains mineral oil, we recommend the addition of 500 µl of Binding Buffer independent of the starting volume. It is also possible to wash the bound PCR fragment once with 500 µl of Binding Buffer.*
- 3. To increase the final DNA yield we recommend an extended incubation time with Elution Buffer (up to 5 minutes), which will lead to a slightly higher final yield.*
- 4. For concentration of PCR-fragments it is possible to elute with lower volume of Elution Buffer than the volume of the starting PCR-mixture. The minimum volume is 10 µl.*

Protocol 4: Purification and concentration of PCR-products over 100 bp with contaminations of primer dimer or shorter fragments

Please read the instructions carefully and conduct the prepared procedure !

for MSB® Vario CleanUp

Note: *Before starting with the purification procedure please place a Spin Filter into a 2.0 ml Receiver Tube !*

1. Binding of the PCR-fragments

Add 10 volumes of **Binding Buffer 2** to 1 volume of the PCR reaction (max. 50 µl) and mix very well by pipetting or vortexing.

Transfer the sample completely onto a Spin Filter and centrifuge for 4 min at maximum speed.

2. Elution of the PCR-fragments

Place the Spin Filter into a new 1.5 ml Receiver Tube.

Add at least the Elution Buffer (or ddH₂O) in **the same volume** as the original PCR reaction (**min. 10 µl**) directly onto the center of the Spin Filter.

Incubate for 1 minute at room temperature. Centrifuge for 1 minute at 10.000 rpm.

Important Notes:

*1. If the PCR-mixture contains mineral oil, we recommend the addition of 500 µl of **Binding Buffer 2** independent of the starting volume. It is also possible to wash the bound PCR-fragment once with 500 µl of Binding Buffer 2.*

*2. To increase the final DNA yield we recommend an extended incubation time with **Elution Buffer** (up to 5 minutes), which will lead to a slightly higher final yield.*

*3. For concentration of PCR-fragments it is possible to elute with lower volume of **Elution Buffer**, than the volume of the starting PCR-mixture. The minimum volume is 10 µl.*

4. At very short PCR-fragments you may get reduced yield so we recommend at least 100 bp in length. Comparing a primer dimer of 70 bp and a PCR-Product of 80 bp there is no big difference. For the purification of short products without contaminations of dimer we recommend the use of protocol 1.

**Especially in DNA-Sequencing and PCR product cloning primer dimer or contamination with short DNA fragments creates problems. T-Vector Cloning kits preferably accept shorter fragments, and the contamination leads to a high cloning background of false clones. In sequencing false fragments like primer dimer may lead to high background in the starting region of the sequence.*

Protocol 5 : Extraction of a DNA-fragment from an agarose gel slice

Please read the instructions carefully and conduct the prepared procedure !

for Invisorb® Spin DNA Extraction Kit and Invisorb® Fragment CleanUp

Important: *TBE-gels contain more potentially inhibitors for down stream application than TAE-gels. So we recommend the use of TAE-gels for critical downstream application!
Before starting with the purification procedure please place a Spin Filter into a 2.0 ml Receiver Tube !*

1. Excise the DNA-fragment from the agarose gel with a sharp scalpel. Minimize the agarose gel slice. Check the weight.
For gel slices up to 150 mg add 500 µl Gel Solubilizer S.
For gel slices > 150 mg add 1 ml of Gel Solubilizer S.
Do not use more than 300 mg gel slice for one Spin Filter.
Transfer the gel slice into a 1.5 or 2.0 reaction tube.
2. Incubate at 50°C for 10 minutes until the agarose gel slice is completely dissolved.
Incubation under continuous shaking (e.g. Eppendorf Thermo mixer) is very helpful.
3. Add 250 µl Binding Enhancer to a 500 µl reaction volume or 500 µl Binding Enhancer to a 1 ml reaction volume and mix the suspension by pipetting some times or by vortexing.
Load appr. 800 µl of the sample onto the Spin Filter. Centrifuge at 10.000 – 12.000 x g for 1 minute. Discard the filtrate. For reaction volumes > 800 µl reload the residual volume onto the Spin Filter and repeat the centrifugation step.
4. Add 500 µl Wash Buffer to the Spin Filter and centrifuge for 30 sec at 8.000 x g
Discard the filtrate. Repeat the washing step once again.
5. Discard the filtrate. Remove the residual ethanol of the Wash Buffer by centrifugation for 4 min at maximum speed (12.000 – 14.000 x g).
6. Transfer the Spin Filter into a new 1.5 ml Receiver Tube.
Add at least 20 µl Elution Buffer directly onto the center of the Spin Filter.
Incubate at room temperature for 5 minutes. Centrifuge for 1 minute at 8.000 x g.

Note: *To increase the final DNA yield we recommend to use a higher volume of Elution Buffer. Please take into account that an increasing volume of Elution Buffer reduces the final concentration of the purified DNA.
An extended incubation time with Elution Buffer (up to 10 minutes) leads also to a slightly higher final yield.*

Protocol 6 : Instruction for repair of A-overhangs in DNA fragments after gel purification

for Invisorb® Spin DNA Extraction Kit and Invisorb® Fragment CleanUp

Transfer 30 µl of the extracted PCR product into a 1.5ml reaction tube.

Add 2 Units Standard Taq DNA Polymerase (no proofreading activity, article number 302030XX), 3.5 µl Taq Buffer, 0,6 µl dNTP's (10 mM each) MgCl₂.
The final concentration in the mixture must be 1.5 mM

Incubate for 15 min at 72°C under continuous shaking in a thermo mixer.

After this treatment, the repaired fragment can be used in cloning/ligation experiments

Troubleshooting of the MSB[®] Spin PCRapace, Invisorb[®] Fragment CleanUp and MSB[®] Vario Cleanup

Problem	Cause	Comments and suggestions
low recovery	poor elution of DNA	add the elution buffer directly onto the centre of the Spin Filter (even if a small elution volume is used).
	problems with mineral oil	apply the correct centrifugation steps Take a higher volume of Binding Buffer Wash once with Binding Buffer

Troubleshooting of the Invisorb[®] Spin DNA Extraction Kit and Invisorb[®] Fragment CleanUp

Problem	Cause	Comments and suggestions
low recovery	incorrect Wash Buffer or no ethanol added	prepare the Wash Buffer exactly as described in the manual. storage of Wash Buffer with firmly fixed cap.
	poor elution of DNA TBE buffered gels are used ineffective solubilization of the agarose gel slice no Binding Enhancer added	add the Elution Buffer directly onto the centre of the Spin Filter (even if a small elution volume is used). the binding of DNA fragments under TBE buffer condition is slightly reduced for smaller fragments than 500 bp please use TAE agarose gels
problems with downstream application, e.g. ligation	contamination with salt components contamination with agarose traces contamination of the final DNA with ethanol	the gel slice must be completely dissolved add the amount of Binding Enhancer needed to the solubilized suspension. Washing of the Spin Filters as described in the manual prolong the incubation time with Wash Buffer to 5 minutes before centrifugation. wash the Spin Filter one time with Gel Solubilizer S . keep the given centrifugation time, extend it if necessary (test the smell)

Order information

Product	Package Size	Order No
MSB [®] Spin PCRapace	10 purifications	1020220900
MSB [®] Spin PCRapace	50 purifications	1020220200
MSB [®] Spin PCRapace	250 purifications	1020220300
MSB [®] Spin PCRapace	500 purifications	1020220400

single components

Binding Buffer	50 ml	1020222000
Elution Buffer	15 ml	1020224000

MSB [®] Vario Cleanup	10 purifications	1020230900
MSB [®] Vario Cleanup	50 purifications	1020230200
MSB [®] Vario Cleanup	250 purifications	1020230300
MSB [®] Vario Cleanup	500 purifications	1020230400

single components

Binding Buffer 1	50 ml	1020232000
Binding Buffer 2	50 ml	1020232100
Elution Buffer	15 ml	1020234000

Invisorb [®] Spin DNA Extraction Kit	3 purifications	1020110100
Invisorb [®] Spin DNA Extraction Kit	50 purifications	1020110200
Invisorb [®] Spin DNA Extraction Kit	250 purifications	1020110300

single components

Gel Solubilizer	30 ml	1020118000
Binding Enhancer	30 ml	1020112900
Wash Buffer (add 42ml ethanol)	18 ml	1020113000
Elution Buffer	15 ml	1020114000

Invisorb [®] Fragment CleanUp	10 purifications	1020300900
Invisorb [®] Fragment CleanUp	50 purifications	1020300200
Invisorb [®] Fragment CleanUp	250 purifications	1020300300

single components

Binding Buffer	50 ml	1020302000
Binding Enhancer	30 ml	1020302900
Wash Buffer (add 42 ml ethanol)	18 ml	1020303000
Elution Buffer	15 ml	1020304000
Gel Solubilizer S	30 ml	1020308000

related products

Invisorb [®] DNA CleanUp	10 purifications	1020400100
Invisorb [®] DNA CleanUp	50 purifications	1020400200
Invisorb [®] DNA CleanUp	250 purifications	1020400300

-