



## SuperHiFi High Fidelity DNA Polymerase 2x Master Mix

Cat. No.	Size Reactions	SuperHiFi High Fidelity DNA Polymerase 2x Master Mix 4 mM MgCl <sub>2</sub>
<b>ID No.</b>		<b>5500400</b>
<b>AO620701</b>	100	2 x 1.25 ml
<b>AO620703</b>	500	10 x 1.25 ml
<b>AO620706</b>	2500	50 x 1.25 ml
<b>AO620707</b>	5000	25 x 5 ml

### Key Features

- Convenient reaction set-up
- High fidelity – measured up to 50x *Taq* fidelity
- Excellent coverage for amplification of difficult amplicons with low to high GC content
- Long range capability: 8.5 kb for gDNA and ≤ 12.5 kb for λDNA
- Recommended for cloning, mutagenesis and other molecular applications requiring extremely high fidelity

SUPERHIFI High Fidelity DNA Polymerase 2x Master Mix is an all-in-one 2x master mix containing the SUPERHIFI High Fidelity DNA Polymerase, SUPERHIFI Buffer, dNTPs and MgCl<sub>2</sub>. Simply mix SUPERHIFI High Fidelity DNA Polymerase 2x Master Mix with primers, DNA template and water and you are ready to carry out PCR.

SUPERHIFI High Fidelity DNA Polymerase is a thermostable DNA Polymerase with proofreading ability. SUPERHIFI High Fidelity DNA Polymerase exhibits 5'→3' DNA polymerase activity and 3'→5' proofreading exonuclease activity. The latter allows the enzyme to correct mis-incorporated nucleotides. The SUPERHIFI High Fidelity DNA Polymerase 2x Master Mix exhibits robust amplification of targets with low to high GC content, as well as a fidelity\* up to 50x *Taq*.

\*Estimated by NGS technology using the Illumina MiSeq instrument.

### Protocol

Reaction conditions such as incubation times, temperatures and amount of template DNA may vary and must be determined individually. Amplification of templates with high GC content, high secondary structures as well as long range amplification may require more optimization.

Set up reaction mixtures in an area separate from that used for DNA preparation or product analysis. **Work on ice at all times.**

1. Thaw SuperHiFi High Fidelity DNA Polymerase 2x Master Mix and primer solutions.  
**It is recommended to completely thaw and thoroughly mix the master mix to ensure proper resuspension of precipitates.**

2. Prepare the reaction mix. Table 1 shows the reaction set up for a final volume of 50 µl. If desired, the reaction size may be scaled down.

**Table 1. Reaction components**

Component	Vol./reaction*	Final concentration*
2x Master Mix	25 µl	1x
Primer A (10 µM)	1 µl	0.2 µM
Primer B (10 µM)	1 µl	0.2 µM
25 mM MgCl <sub>2</sub>	0 µl (0 – 6 µl)	2 mM (2 – 4.5 mM)
Betaine (5M) (optional)	10 – 20 µl	1 – 2M
PCR-grade H <sub>2</sub> O	X µl	-
Template DNA	X µl	genomic DNA: 50 ng (10 – 500 ng) plasmid DNA: 0.5 ng (0.1 – 1 ng) bacterial DNA: 5 ng (1 – 10 ng)
<b>TOTAL volume</b>	<b>50 µl</b>	<b>-</b>

\* Suggested starting conditions; theoretically used conditions in brackets.

3. Mix the reaction mix thoroughly and dispense appropriate volumes into reaction tubes. Mix gently, e.g. by pipetting the master mix up and down a few times.
4. Add template DNA to the individual tubes containing the reaction mix.
5. Program the thermal cycler according to the manufacturer's instructions. For maximum yield and specificity, temperatures and cycling times should be optimized for each new template target or primer pair.
6. Place the tubes in the thermal cycler and start the reaction.

**Table 2. Three-step PCR program**

Cycles	Duration of cycle	Temperature
1	2 – 5 minutes <sup>a</sup>	95 °C
25 – 35	10 – 30 seconds <sup>b</sup> 20 – 40 seconds <sup>c</sup> 30 seconds <sup>d</sup>	95 °C 55 – 70 °C 72 °C
1	5 minutes <sup>e</sup>	72 °C

<sup>a</sup>. Initial denaturation step (optional).

<sup>b</sup>. Denaturation step: This step is the first regular cycling event and consists of heating the reaction to 95 °C for 10 – 30 seconds. It causes melting of the DNA template by disrupting the hydrogen bonds between complementary bases, yielding single-stranded DNA molecules.

<sup>c</sup>. Annealing step: The reaction temperature is lowered to 55 – 70 °C for 20 – 40 seconds allowing annealing of the primers to the single-stranded DNA template. Typically, the annealing temperature is about 3 – 5 °C below the T<sub>m</sub> (melting temperature) of the primers used. **Because of the high salt content within the SUPERHIFI High Fidelity DNA Polymerase 2x Master Mix, annealing temperature will likely be higher than with more traditional PCR master mixes.**

<sup>d</sup>. Extension/elongation step: SuperHiFi High Fidelity DNA polymerase has its optimum activity temperature at 72°C. At this step the DNA polymerase synthesizes a new DNA strand complementary to the DNA template strand. The extension time depends on the length of the DNA fragment to be amplified. **Generally, we recommend an extension time of 1 minute per kb, especially for longer amplicons.**

<sup>e</sup>. Final elongation: This single step is occasionally performed at a temperature of 72 °C for 2 - 5 minutes after the last PCR cycle to ensure that any remaining single-stranded DNA is fully extended.

### Notes for optimization of PCR conditions:

- The optimal MgCl<sub>2</sub> concentration should be determined empirically, but in most cases a final concentration of 2 mM, as provided in SuperHiFi High Fidelity DNA Polymerase 2x Master Mix, will produce satisfactory results. Table 3 provides the volume of 25 mM MgCl<sub>2</sub> to be added to the master mix if a higher MgCl<sub>2</sub> concentration is required.

**Table 3. Additional volume (µl) of MgCl<sub>2</sub> per 50 µl reaction**

Final MgCl <sub>2</sub> conc. in reaction (mM)	2.0	2.5	3.0	3.5	4.0	4.5	5
Volume of 25 mM MgCl <sub>2</sub>	0	1	2	3	4	5	6

- For difficult amplicons, such as GC-rich templates, those with secondary structures or very long amplicons the addition of 1 - 2 M Betaine solution might improve reaction performance (See Additional Products for ordering information). 5 - 10 % DMSO can also be used to improve reaction performance. If using high concentrations of DMSO the annealing temperature has to be lowered as it decreases the primer T<sub>m</sub>.
- Primers of 20 – 40 nucleotides with a GC content of 40 - 60 % are recommended. Online Software such as the Primer3plus <https://primer3plus.com/cgi-bin/dev/primer3plus.cgi> can be used to design primers.

### Kit Components

#### SuperHiFi High Fidelity DNA Polymerase 2x Master Mix

- SuperHiFi High Fidelity DNA Polymerase
- Optimized buffer components, 4.0 mM MgCl<sub>2</sub>
- dNTs

#### 5 M Betaine Enhancer Solution

Sold separately.

Cat No.: AO351104

#### Recommended Storage and Stability

Long term storage at -20 °C. Product expiry at -20 °C is stated on the label.

Optional: Store at +4 °C for up to 6 months.

#### Quality Control

SUPERHIFI High Fidelity DNA Polymerase is tested for contaminating activities with no traces of endonuclease activity or nicking activity.

#### Unit Definition

One unit is defined as the amount of enzyme that incorporates 10 nmoles of dNTPs into an acid-precipitable form of DNA in 30 minutes at 72 °C under standard assay conditions.

### Related Products

SUPERHIFI High Fidelity DNA Polymerase 2x Master Mix	Cat. No.
100 reactions	AO620701
500 reactions	AO620703
2500 reactions	AO620706
5000 reactions	AO620707

SUPERHIFI High Fidelity DNA Polymerase 2 U/µl	Cat. No.
With 10x SUPERHIFI Buffer – 100 U	AO610401
With 10x SUPERHIFI Buffer – 500 U	AO610403
With 10x SUPERHIFI Buffer – 1000 U	AO610404
With 10x SUPERHIFI Buffer – 2500 U	AO610406

PCR Grade Water	Cat. No.
6 x 5 ml	AO360056

Betaine Enhancer Solution 5 M	Cat. No.
5 x 1 ml	AO351104

Reagents for *in vitro* laboratory use only.

Other product sizes, combinations and customized solutions are available. Please look at [www.asone-int.com](http://www.asone-int.com) or ask for our complete product list for PCR Enzymes. For customized solutions please contact us.

Issued 09/2019